

Infusion Pump Analyzer

IPA 2000

New



IPA 2000

IPA 2000

IPA 2000 – Infusion Pump Analyzer is the latest addition to our line of high quality Biomedical test instruments. It is feature rich and ‘affordable’.

IPA 2000 – Infusion Pump Analyzer utilizes our Patented (7669471) technology to quickly and accurately measure flow, volume, backpressure and occlusion pressure without elaborate setup procedures.

IPA 2200 is compact, light weight and easy to use. It is ideal for testing all types of Infusion Pumps. Test results can be easily stored internally and uploaded to a PC via serial port for record keeping.

Features

- Patented Technology for fast and accurate Flow measurement.
- Measures Flow Rate, Volume, Occlusion and Backpressure
- Tests all Infusion pumps including Syringe Pumps
- Test results are automatically saved to internal flash.
- Stored test data can be uploaded to a PC via serial port.
- Large graphics LCD with full function keypad.
- Compact, Lightweight, and Rugged
- Menu Driven and easy to operate
- Two year warranty
- Best value.



Made in USA

Specifications

- Flow Rate: 1 mL - 999 mL/h
- Accuracy: +/- 1%
- Volume: 0 - 999 mL
- Accuracy: +/- 2%
- Occlusion Pressure: 0 – 30 PSI.
- Back Pressure: 0 – 300 mmHg.
- Accuracy: +/- 1%
- Power: 9 V DC Adapter
- Dimension: 6 x 4 x 3 Inches
- Weight: 1 Lbs.



深圳为尔康科技有限公司 联系人：曾祥满 手机：13632925349

QQ：274798107 电话：0755-28896837 地址：深圳市龙岗区沙平北路111号608A

网址：www.medicalQC.com 邮箱：szchina1718@163.com

